Measurements with HyPix-3000, 2D



Seamless Configuration Switching from 2D to 1D, 0D mode





Optics for 2D measurement with HyPix-3000

- PB selection slit + CBO-f + Standard slit box with DHL 0.5 mm
- PB selection slit + CBO-f + PB collimator adapter + Long type Collimator
- PB 0.* mm +PB Collimator adapter + Long/Short type Collimator
- MA + PB Collimator adapter + Long/Short type Collimator
- MA + Standard slit box with DHL 0.5 mm



How to measure the "Wide range Reciprocal Space Map" with HyPix-3000 2D mode



Prepared Example of Macro in SmartLab Guidance software.

Just run from top to bottom like other packages.



Optics Alignment for 2D measurement (Point beam, CBO-f)

Automatic detector distance and center position alignment

Sample height alignment

Wide range RSM measurement

Ex	ec.	Scan axis		Mode	Range	Start (deg)	Stop (deg)	
			▼					
1	1	2-Theta/Omega	-	Continuous 🔻	Absolute 🔻	20.0000	80.0000	Chi :0deg
2	1	2-Theta/Omega	-	Continuous 🔻	Absolute 🔻	20.0000	80.0000	Chi :15deg
3	1	2-Theta/Omega	-	Continuous 🔻	Absolute 🔻	20.0000	80.0000	Chi :30deg
4	1	2-Theta/Omega	-	Continuous 🔻	Absolute 🔻	20.0000	80.0000	Chi :45deg
5		2 Thats Omage	_	Continuous -	Abaaluta –	20.0000	80.0000	

Thin Films Evaluation with HyPix-3000

Epitaxial film : Single crystal film grown on single crystal substrate

 -2D diffraction image analysis

 Orientation information and Mosaic spread (Tilt distribution)

 -Wide range reciprocal space map

 Orientation relations of each layer,

 (and Domains information, unknown phase determination).

 ---- Switch from 2D to 0D configuration

 ---- Twist distribution

 -X-ray reflectivity



Thickness evaluation

In-Plane XRD / Epitaxial film

- Optics alignment
- Sample alignment
 - $\ Z, \, \omega, \, Rx, \, Ry$
- Axis Alignment (InPlane : reflection)
 Set Conditions *use Manual Control dialog box
 - IS : 0.1 mm
 - RS1, RS2 : 20 mm (open)
 - $\omega = 2\theta$: 0.5deg (incident angle)
 - Move to : diffraction angle of investigated material
 - φ scan
 - $\phi: \pm 5 \text{ deg}$
 - Step : 0.1 deg
 - Move to the peak top
- Data Measurement
 - 2θχ/φ scan
 - Φ scan







Divergence angle of incident X-ray beam

Name	Function	Monochromaticity, $\Delta \lambda / \lambda$	Divergence angle, $\Delta \theta$
Multilayer mirror	Suppression of divergence angle	Cu Kα/Kβ intensity ratio from 100:25 to 1000:1	0.04°
Channel-cut monochromator	Suppression of divergence angle Improvement in monochromaticity	3.8 × 10 ⁻⁴	0.0083° (Near $\theta = 0^{\circ}$)
Four-crystal monochromator	Suppression of divergence angle Improvement in monochromaticity	Ge (220) four-crystal: 1.5 $\times 10^{-4}$ Ge(440) four-crystal: 2.3×10^{-5}	Ge (220) four-crystal: 0.0034° Ge (440) four-crystal: 0.0015°



Receiving angle resolution

Name	Function	Slit width (DS, SS, RS)	2θ resolution, $\Delta 2\theta$
		Open	0.84°
		2 mm	0.34°
	Improvement in 2θ resolution	1 mm	0.17°
Slit		0.5 mm	0.084°
		0.2 mm	0.034°
		0.1 mm	0.017°
		0.05 mm	0.0084°
Analyzer crystal	Improvement in 2θ resolution		0.003°
	1		(Often used for thin-
	Improvement in 2θ or		film samples)
Soller slit	$2\theta\chi$ resolution		1°
	(Depends on direction		0.4°
	of insertion)		0.1°



X-rays are useful in daily life

- Soil
- Cosmetics
- Medicine
- Nut
- Plastic bottle
- Fiber
- Juice can
- Pigment









What is a crystal?

- A crystal is defined as a solid composed of atoms arranged in a periodic pattern in three dimensions
- The smallest repeating unit is called a unit cell.

Primitive unit cell





A unit

cell





From the symmetry classification, there are seven crystal systems. The acquired shapes of lattices are classified into 14 types of Bravais lattices.



Crystal Systems	Crystal Axes	Bravais Lattices&Symbol Lette	r Characteristic Patterns
cubic	a=b=c $\alpha=\beta=\gamma=90^{\circ}$	primitive P body-centered I face-centered F	
tetragonal	a=b≠c α=β=γ=90°	primitive P body-centered I	
orthorhombic	a≠b≠c α=β=γ=90°	primitive P body-centered I base-centered C	
rhombohedra trigonal or	l a=b=c α=β=v≠90°	face-centered F primitive P	
hexagonal	$a=b\neq c$ $\alpha=\beta=90^{\circ} \gamma=$	primitive P 120°	
monoclinic	a≠b≠c α=β=90° ≠γ	primitive P base-centered C	
triclinic	a≠b≠c α≠β≠γ ≠90°	primitive P	
2 Rige	aku		۲ 226

Miller indices

- The lattice plane is defined by the three integer numbers (h k l), which are called Miller indices
- The distance between two lattice planes is called lattice spacing





Examples of the Miller indices

Lattice planes and lattice spacing





Examples of Miller indices

 All lattice points can be positioned on flat plane groups denoted by Miller indices





Lattice planes and interplanar spacing

- The **lattice plane** groups of the indicates (*h k l*) are arranged in equal intervals
- This interval is called the interplanar spacing and is denoted as dнк

For the crystal with cubic symmetry:



